

SPECIAL PRICE PROMOTION

## ErgoOne® Manual Pipettes

ErgoOne



### Setting the new standard in pipetting!

STARLAB's variable volume pipettes are designed in accordance with the latest ergonomic principles. ErgoOne® is extremely easy to use and prevents hand fatigue with its minimal operating forces. ErgoOne® provides accurate and precise pipetting in a variety of volume ranges and is available in both single and multi-channel models.


### Features

- Minimal operating forces
- Lightweight ergonomic design
- Colour coded
- Volume lock
- Fully autoclavable
- Low maintenance
- 3-year warranty



**Special Offers!**

See over page for details



[www.starlab.ch](http://www.starlab.ch)

STARLAB Switzerland, Chilefeldstr. 17, CH-5634 Merenschwand, Switzerland  
Tel. +41 (0) 56 664 5982 Fax +41 (0) 56 664 5988 info@starlab.ch www.starlab.ch



\*ErgoOne® is tested in accordance with EN ISO 8655, using STARLAB's TipOne® tips

## SPECIAL PRICE PROMOTION

### Discounted prices on sets of three single-channel pipettes...

ErgoOne® Single-Channel Pipettes							
Cat. No.	Volume Range (Colour Coding)	Volume $\mu$ l	Systematic Error (Inaccuracy) $\pm\%$		Random Error (Imprecision) $\pm\%$		List Price
			8655-2	ErgoOne*	8655-2	ErgoOne*	
S7100-0125	0.1 –2.5 $\mu$ l (Red)	0.1	125	24	75	12	sFr. 285-
		1.25	10	2.5	6.0	1.5	
		2.5	5.0	1.4	3.0	0.7	
S7100-0510	0.5 –10 $\mu$ l (Red)	1.0	12	2.5	8.0	1.8	sFr. 285-
		5	2.4	1.5	1.6	0.8	
		10	1.2	1.0	0.8	0.4	
S7100-0221 (micro cone)	2 –20 $\mu$ l (Red)	2	10	5.0	5.0	1.5	sFr. 285-
		10	2.0	1.2	1.0	0.6	
		20	1.0	1.0	0.5	0.3	
S7100-0220 (standard cone)	2 –20 $\mu$ l (Yellow)	2	10	5.0	5.0	1.5	sFr. 285-
		10	2.0	1.2	1.0	0.6	
		20	1.0	1.0	0.5	0.3	
S7100-1100	10 –100 $\mu$ l (Yellow)	10	8.0	3.0	3.0	1.0	sFr. 285-
		50	1.6	1.0	0.6	0.3	
		100	0.8	0.8	0.3	0.2	
S7100-2200	20 –200 $\mu$ l (Yellow)	20	8.0	2.5	3.0	0.7	sFr. 285-
		100	1.6	1.0	0.6	0.3	
		200	0.8	0.6	0.3	0.2	
S7110-1000	100 –1000 $\mu$ l (Blue)	100	8.0	3.0	3.0	0.6	sFr. 285-
		500	1.6	1.0	0.6	0.2	
		1000	0.8	0.6	0.3	0.2	
S7150-5000	500 –5000 $\mu$ l (Violet)	500	8.0	2.4	3.0	0.6	sFr. 305-
		2500	1.6	1.2	0.6	0.25	
		5000	0.8	0.6	0.3	0.15	

### 20% OFF multi-channel pipettes...

ErgoOne® Multi-Channel Pipettes							
Cat. No.	Volume Range (Colour Coding)	Volume $\mu$ l	Systematic Error (Inaccuracy) $\pm\%$		Random Error (Imprecision) $\pm\%$		List Price
			8655-2	ErgoOne*	8655-2	ErgoOne*	
<b>8-Channel Pipettes</b>							
S7108-0510	0.5 –10 $\mu$ l (Red)	1.0	24	8.0	16	5.0	sFr. 790-
		5.0	4.8	4.0	3.2	2.0	
		10	2.4	2.0	1.6	1.0	
S7108-1100	10 –100 $\mu$ l (Yellow)	10	32	3.0	12	2.0	sFr. 790-
		50	16	1.0	6.0	0.8	
		100	3.2	0.8	1.2	0.3	
S7108-3300	30 –300 $\mu$ l (Green)	30	26.7	3.0	10	1.0	sFr. 790-
		150	5.3	1.0	2.0	0.5	
		300	2.7	0.6	1.0	0.3	
<b>12-Channel Pipettes</b>							
S7112-0510	0.5 –10 $\mu$ l (Red)	1.0	24	8.0	16	5.0	sFr. 890-
		5.0	4.8	4.0	3.2	2.0	
		10	2.4	2.0	1.6	1.0	
S7112-1100	10 –100 $\mu$ l (Yellow)	10	32	3.0	12	2.0	sFr. 890-
		50	16	1.0	6.0	0.8	
		100	3.2	0.8	1.2	0.3	
S7112-3300	30 –300 $\mu$ l (Green)	30	26.7	3.0	10	1.0	sFr. 890-
		150	5.3	1.0	2.0	0.5	
		300	2.7	0.6	1.0	0.3	

Prices exclude VAT. Offer Prices valid until 30.06.12.

Please quote **Offer Ref: 3006ERGO** when ordering

**STARLAB Switzerland**, Chilefeldstr. 17, CH-5634 Merenschwand, Switzerland  
Tel. +41 (0) 56 664 5982 Fax +41 (0) 56 664 5988 info@starlab.ch www.starlab.ch

Buy **ONE SET** of **THREE** single-channels, (any three models), for **sFr. 675.00.** That's only...

**sFr. 225-**  
per pipette!

OR

Buy **TWO** or **MORE SETS** of **THREE** single-channels, (any models) at only...

**sFr. 199-**  
per pipette!

Buy any **8-channel** pipette at only

**sFr. 632-**  
per pipette!

Buy any **12-channel** pipette at only

**sFr. 713-**  
per pipette!

www.starlab.ch

