

# ErgoOne™ Pipettes

New!



## Ergonomic pipetting with ErgoOne


STARLAB's new variable volume pipettes are designed in accordance with the latest ergonomic principles. ErgoOne is extremely easy to use and prevents hand fatigue with its minimal operating forces. ErgoOne provides accurate and precise pipetting in a variety of volume ranges and is available in both single and multi-channel models.

### Features:

- Easy handling with ergonomic design
- Minimal operating forces
- Extremely accurate and precise
- Lightweight
- Low-maintenance high-quality product
- Volume lock
- Fully autoclavable
- Choice of standard or micro cone on the 20µl model
- CE marked for IVD

15%  
OFF

Introductory  
Offer



ErgoOne

Setting the new standard in pipetting

[www.starlab.ch](http://www.starlab.ch)



# ErgoOne™ Pipettes



## Ordering information / Technical Specifications

ErgoOne Pipettes						
Pipette Volume	Cat. No	Colour Coding	Measured Volume	Systematic measurement deviation (inaccuracy)	Coincidental measurement deviation (imprecision)	List Price <b>OFFER PRICE</b>
<b>Single-Channel Pipettes</b>						
0.1-2.5µl	S7100-0125	Red	0.1µl	± 24%	± 12%	sFr. 285- <b>sFr. 242.25</b>
			1.25µl	± 2.5%	± 1.5%	
			2.5µl	± 1.4%	± 0.7%	
0.5-10µl	S7100-0510	Red	1.0µl	± 2.5%	± 1.8%	sFr. 285- <b>sFr. 242.25</b>
			5.0µl	± 1.5%	± 0.8%	
			10µl	± 1.0%	± 0.4%	
2-20µl (micro)	S7100-0221	Red	2.0µl	± 5.0%	± 1.5%	sFr. 285- <b>sFr. 242.25</b>
			10µl	± 1.2%	± 0.6%	
			20µl	± 1.0%	± 0.3%	
2-20µl (standard)	S7100-0220	Yellow	2.0µl	± 5.0%	± 1.5%	sFr. 285- <b>sFr. 242.25</b>
			10µl	± 1.2%	± 0.6%	
			20µl	± 1.0%	± 0.3%	
10-100µl	S7100-1100	Yellow	10µl	± 3.0%	± 1.0%	sFr. 285- <b>sFr. 242.25</b>
			50µl	± 1.0%	± 0.3%	
			100µl	± 0.8%	± 0.2%	
20-200µl	S7100-2200	Yellow	20µl	± 2.5%	± 0.7%	sFr. 285- <b>sFr. 242.25</b>
			100µl	± 1.0%	± 0.3%	
			200µl	± 0.6%	± 0.2%	
100-1000µl	S7110-1000	Blue	100µl	± 3.0%	± 0.6%	sFr. 285- <b>sFr. 242.25</b>
			500µl	± 1.0%	± 0.2%	
			1000µl	± 0.6%	± 0.2%	
500-5000µl	S7150-5000	Violet	500µl	± 2.4%	± 0.6%	sFr. 305- <b>sFr. 259.25</b>
			2500µl	± 1.2%	± 0.25%	
			5000µl	± 0.6%	± 0.15%	
<b>8-Channel Pipettes</b>						
0.5-10µl	S7108-0510	Red	1.0µl	± 8.0%	± 5.0%	sFr. 790- <b>sFr. 671.50</b>
			5.0µl	± 4.0%	± 2.0%	
			10µl	± 2.0%	± 1.0%	
10-100µl	S7108-1100	Yellow	10µl	± 3.0%	± 2.0%	sFr. 790- <b>sFr. 671.50</b>
			50µl	± 1.0%	± 0.8%	
			100µl	± 0.8%	± 0.3%	
30-300µl	S7108-3300	Green	30µl	± 3.0%	± 1.0%	sFr. 790- <b>sFr. 671.50</b>
			150µl	± 1.0%	± 0.5%	
			300µl	± 0.6%	± 0.3%	
<b>12-Channel Pipettes</b>						
0.5-10µl	S7112-0510	Red	1.0µl	± 8.0%	± 5.0%	sFr. 890- <b>sFr. 756.50</b>
			5.0µl	± 4.0%	± 2.0%	
			10µl	± 2.0%	± 1.0%	
10-100µl	S7112-1100	Yellow	10µl	± 3.0%	± 2.0%	sFr. 890- <b>sFr. 756.50</b>
			50µl	± 1.0%	± 0.8%	
			100µl	± 0.8%	± 0.3%	
30-300µl	S7112-3300	Green	30µl	± 3.0%	± 1.0%	sFr. 890- <b>sFr. 756.50</b>
			150µl	± 1.0%	± 0.5%	
			300µl	± 0.6%	± 0.3%	



**DEMO MODELS**  
**TRY**  
before you  
**BUY**  
available on request

### Carousel Stand



Holds up to 8 ErgoOne pipettes.  
Universal! Also holds many other manufacturers' pipettes.

**S7200-0000**

List: s.FR. 175.00

**OFFER s.FR. 148.75**

Tested in accordance with EN ISO 8655, using STARLAB's TipOne tips  
Prices excl. VAT. Offer Prices valid until 30.09.09

### STARLAB Switzerland

Chiliefeldstrasse 17  
CH-5634 Merenschwand  
Switzerland

T: +41 (0)56 664 59 80  
F: +41 (0)56 664 59 88  
info@starlab.ch

STARLAB Switzerland  
is a division of Life Systems Design AG

[www.starlab.ch](http://www.starlab.ch)

