



# Xmas Offers



0.2ml PCR Tubes				
Cat. No	Description	Pack Size	List Price	Offer Price
I1402-8100	0.2ml PCR Tubes, Flat Cap, Natural	1000	sFr. 93.5	<b>sFr. 42.6</b>
I1402-4300	0.2ml PCR Tubes, Domed Cap, Natural	1000	sFr. 104.5	<b>sFr. 42.6</b>




0.5ml PCR Tubes				
Cat. No	Description	Pack Size	List Price	Offer Price
I1405-8100	0.5ml PCR Tubes, Flat Cap, Natural	1000	sFr. 78.3	<b>sFr. 32.3</b>

0.2ml PCR Strips				
Cat. No	Description	Pack Size	List Price	Offer Price
8-Strip Tubes				
I1402-3500	0.2ml 8-Strip PCR Tubes, Natural	125 Strips	sFr. 104.5	<b>sFr. 74.8</b>
8-Strip Caps				
I1400-0900	8-Strip Caps, Flat (Optically Clear)	125 Strips	sFr. 37.0	<b>sFr. 15.5</b>
I1400-0900	8-Strip Caps, Domed	125 Strips	sFr. 37.0	<b>sFr. 15.5</b>



0.2ml PCR Strip with Individual Attached Caps				
Cat. No	Description	Pack Size	List Price	Offer Price
I1402-2900	0.2ml 8-Strip PCR Tubes, Individually Attached Domed Caps	125 Strips	sFr. 122.2	<b>sFr. 84.5</b>

All PCR tubes certified RNase, DNase and DNA free




96-Well PCR Plate				
Cat. No	Description	Pack Size	List Price	Offer Price
4ti-0750	96-Well PCR Plate, Non-Skirted, Black Grid	25	sFr. 90.8	<b>sFr. 81.7</b>
4ti-0751	Strips of 8 flat optical caps (for 4ti-0750)	300 Strips	sFr. 94.9	<b>sFr. 85.4</b>
02-2796-3000	Adhesive PCR Seal	25	sFr. 34.4	<b>sFr. 23.5</b>



Benchtop Coolers				
Cat. No	Description	Pack Size	List Price	Offer Price
E2396-2004	StarChill PCR, Pink-Purple	2	sFr. 147.9	<b>sFr. 95.0</b>
E2396-2006	StarChill PCR, Yellow-Green	2	sFr. 147.9	<b>sFr. 95.0</b>

High performance chiller rack that changes colour to give a visual indication of its temperature. Can maintain a temperature of 4°C for up to 3 hours! Holds 96 x 0.2ml single tubes, strip tubes, or a 96-well PCR plate



PCR Workstation Rack				
Cat. No	Description	Pack Size	List Price	Offer Price
E2396-9649	PCR Workstation with Lid, Mixed Neon	5	sFr. 69.5	<b>sFr. 37.5</b>

Holds 96 x 0.2ml single tubes, strip tubes, or a 96-well PCR plate  
Pack contains one each of neon blue, green, orange, pink, and yellow.

**Offer prices valid until 31 December 2009**



# Xmas Offers - ErgoOne™ Pipettes

**20% OFF**

**ErgoOne**

Setting the new standard in pipetting

## Ordering Information / Technical Specifications

Pipette Volume	Cat. No	Colour Coding	Measured Volume	Systematic measurement deviation (inaccuracy)	Coincidental measurement deviation (imprecision)	List Price	OFFER PRICE
<b>Single-Channel Pipettes</b>							
0.1-2.5µl	S7100-0125	Red	0.1µl	± 24%	± 12%	sFr. 285-	<b>sFr. 228.00</b>
			1.25µl	± 2.5%	± 1.5%		
			2.5µl	± 1.4%	± 0.7%		
0.5-10µl	S7100-0510	Red	1.0µl	± 2.5%	± 1.8%	sFr. 285-	<b>sFr. 228.00</b>
			5.0µl	± 1.5%	± 0.8%		
			10µl	± 1.0%	± 0.4%		
2-20µl (micro)	S7100-0221	Red	2.0µl	± 5.0%	± 1.5%	sFr. 285-	<b>sFr. 228.00</b>
			10µl	± 1.2%	± 0.6%		
			20µl	± 1.0%	± 0.3%		
2-20µl (standard)	S7100-0220	Yellow	2.0µl	± 5.0%	± 1.5%	sFr. 285-	<b>sFr. 228.00</b>
			10µl	± 1.2%	± 0.6%		
			20µl	± 1.0%	± 0.3%		
10-100µl	S7100-1100	Yellow	10µl	± 3.0%	± 1.0%	sFr. 285-	<b>sFr. 228.00</b>
			50µl	± 1.0%	± 0.3%		
			100µl	± 0.8%	± 0.2%		
20-200µl	S7100-2200	Yellow	20µl	± 2.5%	± 0.7%	sFr. 285-	<b>sFr. 228.00</b>
			100µl	± 1.0%	± 0.3%		
			200µl	± 0.6%	± 0.2%		
100-1000µl	S7110-1000	Blue	100µl	± 3.0%	± 0.6%	sFr. 285-	<b>sFr. 228.00</b>
			500µl	± 1.0%	± 0.2%		
			1000µl	± 0.6%	± 0.2%		
500-5000µl	S7150-5000	Violet	500µl	± 2.4%	± 0.6%	sFr. 305-	<b>sFr. 244.00</b>
			2500µl	± 1.2%	± 0.25%		
			5000µl	± 0.6%	± 0.15%		
<b>8-Channel Pipettes</b>							
0.5-10µl	S7108-0510	Red	1.0µl	± 8.0%	± 5.0%	sFr. 790-	<b>sFr. 632.00</b>
			5.0µl	± 4.0%	± 2.0%		
			10µl	± 2.0%	± 1.0%		
10-100µl	S7108-1100	Yellow	10µl	± 3.0%	± 2.0%	sFr. 790-	<b>sFr. 632.00</b>
			50µl	± 1.0%	± 0.8%		
			100µl	± 0.8%	± 0.3%		
30-300µl	S7108-3300	Green	30µl	± 3.0%	± 1.0%	sFr. 790-	<b>sFr. 632.00</b>
			150µl	± 1.0%	± 0.5%		
			300µl	± 0.6%	± 0.3%		
<b>12-Channel Pipettes</b>							
0.5-10µl	S7112-0510	Red	1.0µl	± 8.0%	± 5.0%	sFr. 890-	<b>sFr. 712.00</b>
			5.0µl	± 4.0%	± 2.0%		
			10µl	± 2.0%	± 1.0%		
10-100µl	S7112-1100	Yellow	10µl	± 3.0%	± 2.0%	sFr. 890-	<b>sFr. 712.00</b>
			50µl	± 1.0%	± 0.8%		
			100µl	± 0.8%	± 0.3%		
30-300µl	S7112-3300	Green	30µl	± 3.0%	± 1.0%	sFr. 890-	<b>sFr. 712.00</b>
			150µl	± 1.0%	± 0.5%		
			300µl	± 0.6%	± 0.3%		



### Features:

- Volume display positioned so always visible during pipetting
- Minimal operating forces
- Extremely accurate and precise
- Lightweight
- Low-maintenance high-quality product
- Volume lock
- Fully autoclavable
- Choice of standard or micro cone on the 20µl model
- CE marked for IVD



Tested in accordance with EN ISO 8655, using STARLAB's TipOne tips  
Prices excl. VAT. Offer Prices valid until 31.12.09

**STARLAB Switzerland**  
Chilefeldstrasse 17  
CH-5634 Merenschwand  
Switzerland

T: +41 (0)56 664 59 80  
F: +41 (0)56 664 59 88  
info@starlab.ch

STARLAB Switzerland  
is a division of Life Systems Design AG

[www.starlab.ch](http://www.starlab.ch)

