Stainless steel freezer racking.

IT'S COOL TO BE ORGANISED

MADE FROM

50/0

RECYCLED

MATERIAL







Stainless steel freezer racking.

QUALITY RACKING UNITS TO OPTIMISE YOUR FREEZER CAPACITY

For use with all standard storage boxes up to 136 x 136 x 53 mm

Racking with pull out shelves also available for use with our popular StarStore 100 boxes.



UPRIGHT FREEZER RACKS, SIDE ACCESS

> Side access shelving system for optimising space in your upright freezer



UPRIGHT FREEZER RACKS, PULL OUT SHELVES

> Pull out shelves provide the convenience of access to your samples without having to remove the rack from the freezer



CHEST FREEZER TOWER RACKS

> Tower racks for optimising storage space and easy access to your samples from your chest freezer





BE COOL, GET ORGANISED!

Starlab's new range of freezer racking provides a lightweight but sturdy solution to keep your samples safe. Manufactured to the highest specifications, laser cutting technology results in precision racking with no sharp edges.

Keeping your freezer tidy and optimising storage space is an efficient use of energy. Our racks are made with 85% recycled stainless steel, in a factory powered 100% by sustainable energy.

So get your samples organised with Starlab's new range of freezer racking. To help, we have even stamped the article number on the rack to make reordering more racks easier. Oh, and we have very short lead times too. Now that is organised!



SMOOTH

Laser cutting technology results in high quality units with no sharp edges



BRANDED

Article No. etched onto every rack for easy reordering

SUSTAINABILITY JUST GOT COOLER!

- Freezer racking helps optimise capacity by keeping freezers tidy and maximising storage space for efficient energy use
- Do not store anything on top of your freezer as this can prevent good air circulation
- Made from 85% recycled stainless steel
- Manufactured in a factory powered by 100% sustainable energy
- Reduced carbon footprint





Stainless steel racking for upright freezers.

1 RACKS WITH SIDE ACCESS





PULL OUT SHELVES

RACKS WITH

Capacity	Format	Ext. Dimensions (WxDxH)	Art. No								
9-Place	3 deep x 3 high	139 x 426 x 169 mm	E6000-0033								
12-Place	3 deep x 4 high	139 x 426 x 224 mm	E6000-0034								
15-Place	3 deep x 5 high	139 x 426 x 297 mm	E6000-0035								
12-Place	4 deep x 3 high	139 x 565 x 169 mm	E6000-0043								
16-Place	4 deep x 4 high	139 x 565 x 224 mm	E6000-0044								
20-Place	4 deep x 5 high	139 x 565 x 279 mm	E6000-0045								
24-Place	4 deep x 6 high	139 x 565 x334 mm	E6000-0046								
Racking suitab	Racking suitable for use with storage boxes up to maximum (WxDxH) 136 x 136 x 53 mm.										




For use with all storage boxes up to 136 x 136 x 53 mm.



J.	

Ä	starlab.click /shelfracks

Capacity	Format	Ext. Dimensions (WxDxH)	Art. No							
9-Place	3 deep x 3 high	140 x 450 x 173 mm	E6001-0033							
12-Place	3 deep x 4 high	140 x 450 x 231 mm	E6001-0034							
15-Place	3 deep x 5 high	140 x 450 x 288 mm	E6001-0035							
12-Place	4 deep x 3 high	140 x 569 x 174 mm	E6001-0043							
16-Place	4 deep x 4 high	140 x 569 x 231 mm	E6001-0044							
20-Place	4 deep x 5 high	140 x 569 x 281 mm	E6001-0045							
24-Place	4 deep x 6 high	140 x 569 x 337 mm	E6001-0046							
28-Place	4 deep x 7 high	140 x 569 x 403 mm	E6001-0047							
Racking suitable for use with storage boxes up to maximum (WxDxH) 136 x 136 x 53 mm.										



Storage boxes are sold separately and are not included with any of the freezer racking.

3 RACKS WITH PULL OUT SHELVES

For use with StarStore 100 storage boxes





Capacity	Format	Ext. Dimensions (WxDxH)	Art. No							
9-Place	3 deep x 3 high	147 x 500 x 179 mm	E6001-0133							
12-Place	3 deep x 4 high	147 x 500 x 239 mm	E6001-0134							
16-Place	4 deep x 4 high	147 x 655 x 239 mm	E6001-0144							
Racking suitable for use with StarStore 100 storage boxes (WxDxH) 140 x 155 x 55 mm.										





Racking guide: upright freezers.

			ЕРР	ENDO	RF (NE	W BR	UNSW	ICK)		PHCBI (PANASONIC / SANYO) - MDF RANGE										THERMO											
	RACKING A SIDE ACCE STANDARE STORAGE I up to 136 x 13	SSS FOR D BOXES	F570h CryoCube	F740h CryoCube	U360VIP	U535VIP	U570	U725VIP	DU302VX	DU502VH-PH	DU502VX	DU702VH-PE	DU702VX	DU900V-PE	U33V	U4186S	U500VX	U53V	U5386S	U55V	U700VX	U7386S	U76V	DW86L288	DW86L338(J)	DW86L388(A/J)	DW86L486	DW86L490(J)	DW86L628	DW86L728(J)	ExF60086
	E6000-0033	9-place (3d x 3h)																													
	E6000-0034	12-Place (3d x 4h)														8*															
ACCESS	E6000-0035	15-Place (3d x 5h)														8*								12		16					
AC(E6000-0043	12-Place (4d x 3h)																													
SIDE	E6000-0044	16-Place (4d x 4h)	25				25								6*				16						6*			8*			
	E6000-0045	20-Place (4d x 5h)							12					14*	6*		8*	8*		8*	12*	24	12*		6*		16	8*	20	24	
	E6000-0046	24-Place (4d x 6h)								16	16	24	24	14*			8*	8*		8*	12*		12*								
(2)	E6001-0033	9-place (3d x 3h)																													
	E6001-0034	12-Place (3d x 4h)														8*															
VES	E6001-0035	15-Place (3d x 5h)														8*								12		16					
SHELVES	E6001-0043	12-Place (4d x 3h)																													
OUT	E6001-0044	16-Place (4d x 4h)	25				25								6*				16						6*			8*			
PUL (E6001-0045	20-Place (4d x 5h)							12					14*	6*		8*	8*		8*	12*	24	12*		6*		16	8*	20	24	
- 4	E6001-0046	24-Place (4d x 6h)								16	16	24	24	14*			8*	8*		8*	12*		12*								
	E6001-0047	28-Place (4d x 7h)		18	9	12		18																							

RACKING WITH PULL OUT SHELVES FOR STARSTORE 100 BOXES

140 x 155 x 55 mm




Table shows the maximum number of racks that fit into the model of freezer.



^{*}These freezers hold racks of different sizes on different shelves; table shows maximum number of each rack size. The maximum number of racks per freezer is the total of each added together. Version 1: more info to be added as soon as it becomes available.


Stainless steel racking for chest freezers.





For use with standard storage boxes up to 136 x 136 x 53 mm

Capacity	Format	Ext. Dimensions (WxDxH)	Art. No						
9-Place	9 high	141 x 148 x 501 mm	E6010-0009						
10-Place	10 high	141 x 148 x 556 mm	E6010-0010						
11-Place	11 high	141 x 148 x 611 mm	E6010-0011						
12-Place	12 high	141 x 148 x 666 mm	E6010-0012						
13-Place	13 high	141 x 148 x 721 mm	E6010-0013						
Racking suitable for use with storage boxes up to maximum (WxDxH) 136 x 136 x 53 mm.									



HOW TO MEASURE RACKING FOR CHEST FREEZERS

starlab.click/chestracks

Racking guide: chest freezers.

			EPPEN	DORF (NE	W BRUNS	SWICK)	РНСЕ	BI (PANAS - MDF	ARTICO	THERMO		
	RACKING F STANDARD STORAGE E up to 136 x 130	C340	C585VIP	0990	C760VIP	C2156-VAN	394	594	794	ULTF 420	CXF ULT2090	
4	E6010-0009	9-high						21				
	E6010-0010	10-high					14					
	E6010-0011	11-high									24	
	E6010-0012	12-high										33
	E6010-0013	13-high	13	13	13	13			24	36		
	Table shows th	ne maximum	number o	f racks that	fit into the	model of fi	eezer.				•	



Storage boxes for freezer racking.

FOR USE WITH FREEZER RACKING THAT HOLDS STANDARD BOXES UP TO 136 x 136 x 53 mm



StarStore 81

81-place polypropylene storage box with transparent lift-off lid. Printed reference grid on base. Natural or neon colours.

WxDxH: 130 x 130 x 45 mm



StarStore 49

49-place polypropylene storage box with transparent lift-off lid. Natural colour only.

WxDxH: 130 x 130 x 45 mm

starlab.click/ppstorage



StarStore

Polypropylene storage box with no dividers. Transparent lift-off lid. Natural colour only.

WxDxH: 130 x 130 x 45mm



Cardboard storage

Moisture resistant cardboard boxes with a choice of gid insert configurations. Boxes also available without grids.

WxDxH: 133 x 133 x 53 mm

starlab.click/cardboardstorage



StarRack 64

64-place polypropylene double-sided rack with transparent lift-off lid. Neon colours.

WxDxH: 130 x 130 x 50 mm

starlab.click/starrac



81-place coloured cryo boxes

Polycarbonate storage boxes with transparent lift-off lid. 9x9 fixed, grid dividers. Printed reference grid on base.

WxDxH: 133 x 133 x 49 mm



100-place coloured cryo boxes

Polycarbonate storage boxes with transparent lift-off lid. 10x10 fixed, prong dividers. Printed reference grid on base.

WxDxH: 133 x 133 x 49 mm



81- and 100-place cryo boxes

Polycarbonate storage boxes. Transparent lift-off lid with printed reference grid. 9x9 box, fixed grid dividers; 10x10 fixed pronged dividers. White base.

WxDxH: 130 x 130 x 55 mm



StarStore 81Z

Zylar resin storage box with transparent lift-off lid. 9x9 fixed, grid dividers. Natural colour only.

WxDxH: 130 x 130 x 55 mm



StarStore 100

100-place polypropylene storage box with fixed 10x10 dividers. Hinged, lockable lid. Positioning lugs on lid for safe stacking. Natural or neon colours.

WxDxH: 140 x 155 x 55 mm

For use with the freezer racks that hold StarStore 100 storage boxes only.

starlab.click/starstore100





starlab

Starlab International GmbH

Neuer Höltigbaum 38 22143 Hamburg Germany

T: +49 (0)40 675 99 39 0 F: +49 (0)40 675 99 39 20 info@starlab.de www.starlabgroup.com

Starlab GmbH

Neuer Höltigbaum 38 22143 Hamburg Germany

T: +49 (0)40 675 99 39 0 F: +49 (0)40 675 99 39 20 info@starlab.de www.starlab.de

Starlab (UK), Ltd

5 Tanners Drive Milton Keynes MK14 5BU United Kingdom

T: +44 (0)1908 283800 F: +44 (0)1908 283802 infoline@starlab.co.uk www.starlab.co.uk

Starlab France SARL

30 Rue Jean Rostand 91400 Orsay France

T: +33 (0)1 60 13 71 70 F: +33 (0)1 69 41 48 65 info@starlab.fr www.starlab.fr

Starlab S.r.l.

Via Pinturicchio 1 20133 Milano Italy

T: +39 (0)2 7020 1040 F: +39 (0)2 7020 1033 info@starlab.it www.starlab.it



IT'S COOL TO BE ORGANISED.

Starlab's new range of freezer racking provides a lightweight but sturdy solution to keep your samples safe. Manufactured to the highest specifications, laser cutting technology results in precision racking with no sharp edges.

Keeping your freezer tidy and optimising storage space is an efficient use of energy. Our racks are made with 85% recycled stainless steel, in a factory powered 100% by sustainable energy.

Get organised, with Starlab's new range of freezer racking!

