

**Trusted for consistency,  
designed for growth.**  
**You are passionate  
for science.**

**SPECIAL OFFER HIGHLIGHTS**

Valid until 31.08.2026



SCAN TO  
BUY ONLINE

**IN THIS EDITION:**  
**FLASKS, DISHES & PLATES**  
**BOTTLE TOP FILTRATION**  
**ASPIRATOR**  
**TUBES & VIALS FOR CELL CULTURE**



## CytoOne® Dishes.

FROM:  
CHF 67.25\*



CytoOne® TC-treated dishes feature stacking rings with vented design for thermal equalisation. Untreated lids reduce condensation: ideal for microscopy.

Size	Pack Size	Art. No.	OFFER*
<b>TC-Treated Dishes</b>			
35 x 10 mm	30 x 10	CC7682-3340	CHF 68.42
60 x 15 mm	30 x 10	CC7682-3359	CHF 67.25
100 x 20 mm	30 x 10	CC7682-3394	CHF 112.64
150 x 20 mm	12 x 5	CC7682-3614	CHF 66.43

Certified and sterile (SAL 10<sup>-9</sup>).

All CytoOne® products are certified RNase, DNase, DNA & Pyrogen free and Gamma sterilised (SAL 10<sup>-9</sup>). TC-Treated for optimum cell attachment and growth.

## CytoOne® Plates.

FROM:  
CHF 77.78\*



Complete 360° chimney well design prevents cross-contamination and ensures a uniform temperature transfer across all wells to prevent the edge effect. Integrated vents and stacking rims create space for continuous airflow between plates and reduces evaporation.

Wells	Pack Size	Art. No.	OFFER*
<b>TC-Treated Plates</b>			
6-Well	50	CC7682-7506	CHF 77.78
12-Well	50	CC7682-7512	CHF 87.48
24-Well	50	CC7682-7524	CHF 89.76
48-Well	50	CC7682-7548	CHF 94.78
96-Well	50	CC7682-7596	CHF 89.18

Certified and sterile (SAL 10<sup>-9</sup>).



## CytoOne® Flasks.

FROM:  
CHF 92.87\*



All flasks supplied in convenient 'mini-packs' in easy-open, resealable zip top bags that protect unused flasks.



Vented filter caps with a 0.2 µm hydrophobic membrane for gas exchange and to exclude bacteria and fungal contaminants.



Two-position plug seal caps allow you to choose secure manual venting or an airtight seal.

Size	Pack Size	Art. No.	OFFER*
<b>TC Treated Flasks, Vented Filter Caps</b>			
T-25	30 x 10	CC7682-4825	CHF 296.91
T-75	20 x 5	CC7682-4875	CHF 167.55
T-150	8 x 5	CC7682-4815	CHF 108.74
T-225	5 x 5	CC7682-4822	CHF 110.23

Certified and sterile (SAL 10<sup>-9</sup>).

Size	Pack Size	Art. No.	OFFER*
<b>TC Treated Flasks, Plug Seal Caps</b>			
T-25	30 x 10	CC7682-4325	CHF 254.34
T-75	20 x 5	CC7682-4175	CHF 140.33
T-150	8 x 5	CC7682-4415	CHF 92.87
T-225	5 x 5	CC7682-4225	CHF 99.89

Certified and sterile (SAL 10<sup>-9</sup>).

## Centrifuge Tubes.

Safeguard your cell viability. Designed to support your most sensitive cell work. Graduated to full 15 ml and 50 ml volumes, with extra graduations on the conical base.

Extra moulded graduations!  
15 ml Conical Tube: 0 – 1.0 ml  
50 ml Conical Tube: 0 – 3.0 ml



Description	Pack Size	Art. No.	OFFER*
15 ml Conical, Loose	20 x 25	E1415-0100	CHF 93.28
15 ml Conical, Racked	20 x 25	E1415-0700	CHF 99.02
50 ml Conical, Loose	20 x 25	E1450-0100	CHF 139.31
50 ml Conical, Racked	20 x 25	E1450-0700	CHF 156.79
50 ml Skirted, Loose	20 x 25	E1450-0400	CHF 154.73

Certified and sterile (SAL 10<sup>-9</sup>).

### MICROTUBE PESTLE

for 1.5/2.0 ml tubes.  
Assorted colours.  
Certified. Pack 10.  
I1415-5390  
CHF 6.68\*



### CULTURE TUBES

14 ml tubes with dual-position caps that can be left loose for aerobic work or sealed for anaerobic culture.

Polypropylene tube  
-190 °C / 3000 xg.  
20 x 25-Tube bags  
I1485-0810  
CHF 127.55\*

Polystyrene tube  
-10 °C / 1400 xg  
20 x 25-Tube bags  
I1485-2810  
CHF 134.73\*



## Serological Pipettes.



Vol.	Qty	Art. No.	OFFER*
1.0 ml	500	E4860-0001	CHF 42.50
2.0 ml	500	E4860-0002	CHF 48.56
5. ml	200	E4860-0005	CHF 27.20
10 ml	200	E4860-0010	CHF 32.05
25 ml	150	E4860-0025	CHF 47.36
50 ml	100	E4860-0050	CHF 44.51
100 ml	50	E4860-1100	CHF 59.20

Certified RNase, DNase, DNA, Latex and Pyrogen free. Sterile (SAL 10<sup>-9</sup>). Individually wrapped.



Replaceable rechargeable battery!

### PIPETTE CONTROLLER.

Pressure-responsive speed controls make the ErgoOne® FAST ideal for cell culture work such as resuspension of bacteria or cell pellets, nucleic acid precipitates and aspiration of a cell layer.

### ErgoOne® FAST

Art. No. S7166-0010  
CHF 276.36\*

## Cryogenic Vials.

Supplied in convenient, self-standing resealable pouches. Certified. Sterile (SAL 10<sup>-9</sup>).



	Vol:	1.0 ml	1.8 ml	1.8 ml	3.6 ml	3.6 ml	4.5 ml	4.5 ml
	Base:	Skirted	Conical	Skirted	Conical	Skirted	Conical	Skirted
<b>EXTERNAL Threaded Tube Self-Sealing Cap, Sterile</b>	<b>E3090-6212</b>	<b>E3090-6321</b>	<b>E3090-6222</b>					<b>E3090-6242</b>
	10 x 50	10 x 50	10 x 50	n/a	n/a	n/a	n/a	6 x 50
	CHF 104.57*	CHF 112.68*	CHF 109.37*					CHF 77.85*
<b>INTERNAL Threaded Tube, Silicone Seal in Cap, Sterile</b>	<b>E3110-6112</b>	<b>E3110-6121</b>	<b>E3110-6122</b>	<b>E3110-6131</b>	<b>E3110-6132</b>	<b>E3110-6141</b>	<b>E3110-6142</b>	
	10 x 50	10 x 50	10 x 50	8 x 50	8 x 50	6 x 50	6 x 50	
	CHF 139.59*	CHF 139.64*	CHF 140.17*	CHF 135.10*	CHF 129.63*	CHF 106.80*	CHF 106.80*	

**\*OFFER PRICES REFLECT A 40% DISCOUNT vs LIST PRICE**

Offer valid to 31.08.26 / while stock lasts.

# ErgoOne® Starter Packs.



ErgoOne® pipettes paired with Starlab's trusted TipOne® pipette tips create the perfect pipetting partnership!

3 x ErgoOne® Single-Channel Pipettes: 20 µl, 200 µl, 1000 µl  
3 x racks of sterile TipOne® Tips  
3 x Shelf-Adhesive Pipette Holders  
1 x Laboratory Timer

**STARTERPACK**  
**CHF 540.40**

4 x ErgoOne® Single-Channel Pipettes: 2.5 µl, 20 µl, 200 µl, 1000 µl,  
3 x racks of sterile TipOne® Tips  
4 x Shelf-Adhesive Pipette Holders  
1 x Laboratory Timer

**STARTERPACK42.5**  
**CHF 716.58**

4 x ErgoOne® Single-Channel Pipettes: 10 µl, 20 µl, 200 µl, 1000 µl,  
3 x racks of sterile TipOne® Tips  
4 x Shelf-Adhesive Pipette Holders  
1 x Laboratory Timer

**STARTERPACK410**  
**CHF 716.58**

## Cell Scrapers.

Gently raise/lift the cells with Starlab cell scrapers. The silicone blades of 220 mm and 250 mm scrapers prevent cell lysing. The flexible ABS handles on all scrapers provides easy handling inside flasks. Individually wrapped. Certified. Sterile.



220 mm straight handle, 11 mm fixed blade  
CC7600-0220  
Pack Size: 100  
SAL 10<sup>5</sup> **CHF 212.98\***



250 mm angled handle, 11 mm pivoting blade  
CC7600-0250  
Pack Size: 100  
SAL 10<sup>5</sup> **CHF 157.42\***



300 mm straight handle, 20 mm pivoting blade  
CC7600-0300  
Pack Size: 150  
SAL 10<sup>6</sup> **CHF 147.88\***

## Win a TipOne® lego set!

Enter our online competition for the chance to win!

Scan the QR code or visit the link.

Competition closes: 31.08.2026



[starlab.click/win-lego](http://starlab.click/win-lego)

**\*OFFER PRICES REFLECT A 40% DISCOUNT vs LIST PRICE**




Offer valid to 31.08.26 / while stock lasts.

## Bottle Top Filtration.

**Flexible system!** Choose complete full assembly units for ready-to-use filtration, GL45 compatible bottle top funnels to pair with your own bottles or media bottles for versatile storage and preparation.

Superior filtration performance!



Bottle Top Filtration System					
	150 ml	250 ml	500 ml	1000 ml	
<b>Full Assembly Unit</b>					
	0.22 µm PES membrane Pack size: 12	<b>CC8221-3226</b> <b>CHF 80.27*</b>	<b>CC8221-4226</b> <b>CHF 89.52*</b>	<b>CC8221-5226</b> <b>CHF 104.44*</b>	<b>CC8221-6226</b> <b>CHF 166.90*</b>
<b>Funnel Tops only (GL45 compatible)</b>					
	0.22 µm PES membrane Pack size: 24	<b>CC8222-3226</b> <b>CHF 85.06*</b>	<b>CC8222-4226</b> <b>CHF 95.93*</b>	<b>CC8222-5226</b> <b>CHF 122.19*</b>	<b>CC8222-6226</b> <b>CHF 199.23*</b>
<b>Media Bottles only</b>					
	Bottle and Cap Pack size: 12	<b>CC8221-3000</b> <b>CHF 61.63*</b>	<b>CC8221-4000</b> <b>CHF 69.49*</b>	<b>CC8221-5000</b> <b>CHF 77.08*</b>	<b>CC8221-6000</b> <b>CHF 120.85*</b>

All products are certified RNase, DNase and Pyrogen free, are sterile (SAL 10<sup>-6</sup>), and supplied individually wrapped.

### THE PERFECT PARTNER

Starlab Aspirator. N2400-9000.  
**CHF 1,283.43\***



## Syringe Filters.

Designed for efficiency and safety. Low hold-up volume means you can recover more of your valuable solution. 33 mm diameter filters supplied in easy-open blister packs showing sterile validity date (SAL 10<sup>-6</sup>).

100 Filters: Double the usual pack size!



0.22 µm PES Sterile.  
E4870-1226  
Pack Size: 100  
**CHF 111.47**



0.22 µm PVDF Sterile.  
E4870-1221  
Pack Size: 100  
**CHF 117.34\***



0.22 µm CA Sterile.  
E4870-1223  
Pack Size: 100  
**CHF 111.47**

## Cell Strainers.

Cell Strainers provide a simple and reliable way to achieve uniform cell suspensions. Designed to fit perfectly into 50 ml centrifuge tubes. Gamma sterilised (SAL 10<sup>-6</sup>). Individually wrapped.



40 µm mesh. Sterile.  
CC8111-0042  
Pack Size: 50  
**CHF 78.42\***



70 µm mesh. Sterile.  
CC8111-0072  
Pack Size: 50  
**CHF 78.42\***



100 µm mesh. Sterile.  
CC8111-0102  
Pack Size: 50  
**CHF 78.42\***